












CONTROL BUTTONS



1.  **On - Off**
2.  **Timer.** Waveform cycle starts when electrodes touch skin. Countdown timers: select 20, 40 or 60 minutes.
3.  **Lock.** Press 2 times within 5 seconds to lock or unlock settings during treatment if desired.
4.  or  **Current.** Increases or decreases current. Increases at 50 microamperes (μA) per second when held down. Decreases at 100 μA per second. When the device is turned on, the current defaults to 100 μA .

LCD DISPLAY

- Light Sensor** lights LCD screen for 10 seconds in a dark room when any button is pushed.
-  **Test Circuit** symbol and an audio warning indicates device is not treating. The timer will stop and the device will turn off in 30 minutes if nothing else is done. Current must be set above 0 μA and moistened electrodes must be in contact with skin for Alpha-Stim[®] AID to work. An audible signal and resumption of the timer occurs the instant the integrity check circuitry determines everything is working properly.
-  **Battery Charge** indicator. Replace battery when only 1 bar remains at which time the device shall give out an audible low battery warning and repeat the warning every 10 minutes (unless mute function is on).
-  **Timer**. Select 20, 40 or 60 minutes countdown to auto-off.
- Time** setting indicator. Displays time remaining.
-  0 - 500 microamperes (μA) of current.
-  **Mute** feature. To turn on or off press the key sequence: **Lock-Timer-Lock**.
-  **Lock** is on when display is lit.
- Σ Cumulative Timer** records total hours and minutes Alpha-Stim[®] AID has been in use.

